

Updated 8/1/2014



**Camera Button Utility Chart**

Software	Mode	Time	Key 1	Key 2	Key 3	Notes
Apteryx XVLite or XRayVision	Send All Mode	50ms	SPACE			Includes other OEM versions of XV and XVLite Only needed for QuickCam Duo and HD
Curve Hero	Send All Mode	50ms	ENTER			
Dentrix Image	Send All Mode	50ms	ALT+C	ALT+S		
Dexis	Send All Mode	50ms	T	SPACE		Dexis 10 only requires Key 1 to be set.
EagleSoft	Send All Mode	50ms	F	S		
EasyDental Imaging	Send All Mode	50ms	ALT+C	ALT+S		
Kodak Dental Imaging	Send All Mode	50ms	F5	F5		
Mogo	Send All Mode	50ms	F8			
Schick CDR Dicom	Send All Mode	50ms	ALT+C	U		
SoftDent Imaging	Send All Mode	50ms	ALT-S			
VixWin Platinum	Send All Mode	50ms	SPACE	ALT+G	SPACE	
VixWin Pro	Send All Mode	50ms	ALT+F	ALT+G	ALT+F	
XDR	Send All Mode	50ms	CTRL+I	CTRL+I		